

DIRECT OPHTHALMOSCOPE

YZ6H CE

Direct Ophthalmoscope

- LED cold light source.It is durable enough with the designed service life of 100000 hours.
- Low heat benefits from the cold light source.
- Bulb maintains 3200K color temperature for long time.The illumination is soft and more similar with sunlight.
- Ultraviolet radiation is eliminated to fully protect patient' s eyes.
- Even light spot and sharp observed fundus image.
- Soft rubber around the eyeholes protects eyes and stabilizes the ophthalmoscope.
- Aseismic design,soft start and automatic turn-off function ensure the product to be safer and more durable.



Specifications

Illumination Type	Large Spot , Small Spot , Slit , Net , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	3.5V/1W LED Bulb
Power Supply	AC110V/220V±10% , 50Hz/60Hz

YZ6F CE

Direct Ophthalmoscope

- Excellent optical performance and continuously adjustable brightness, easy to observe the erect fundus image.
- Germany bulb provides the high color temperature,good color rendering,even light spot and sharp image.
- Complete types of light spots and filters meet various observation requirements.
- Soft rubber around the eyeholes protects user' s eyes and stabilizes the ophthalmoscope.
- New functions of soft starting and automatic turn-off extend the service life of bulb.
- Convenient to place and easy to operate.



Specifications

Illumination Type	Large Spot , Small Spot , Slit , Net , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	3.5V/2.8W Halogen Bulb
Power Supply	AC110V/220V±10% , 50Hz/60Hz

DIRECT OPHTHALMOSCOPE

YZ11 FDA CE

Direct Ophthalmoscope

- Elegant appearance.
- It is convenient for mobile inspection.
- Both dry battery and rechargeable battery power are available.



Specifications

Illumination Type	Large Spot , Middle Spot , Small Spot , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	2.6V/0.78W , Incandescent Bulb
Power Consumption	C Alkaline battery , rechargeable battery

YZ11D FDA CE

Rechargeable Direct Ophthalmoscope

- Improved optical design and excellent optical performance.
- Imported halogen bulb provides high and adjustable brightness and good color rendering.
- Smart, portable and convenient for mobile inspection.
- Equipped with special charger quickly and conveniently recharging, against over recharging.
- Anti-recharging and anti-discharging function protects the battery.



Specifications

Illumination Type	Large Spot , Small Spot , Slit , Net , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	3.5V/2.8W Halogen Bulb
Power Supply	AC110V/220V±10% , 50Hz/60Hz
Power (Main Body of Ophthalmoscope)	Li Battery Pack 3.7V , Rechargeable