

OPERATION MICROSCOPE



YZ20P Operation Microscope multi-section

- The whole microscope is smart.
- Three-step magnifications.
- Spring balance system in the arm ensures the microscope move freely.
- All of the lens are multi-coated, mildewproof and anti-reflection.
- For optional objectives and other parts applicable to ENT, surgery and Gynaecology.



Specifications

Magnification of objectives	12.5×	
Adjustable range of diopter	± 6D	
Objective	200mm	300mm(Optional)
Vision field	6.6×(Φ35mm) 10×(Φ22mm) 15×(Φ14mm)	4.2×(Φ52mm) 6.7×(Φ34mm) 11×(Φ20mm)
Adjustment Range of Pupil Distance	50mm ~ 70mm	
Illumination Source	12V/100W ,Halogen Lamp	
Illumination Type	6° Coaxial Illumination	
Reaching Range of Arm	860mm	
Vertical movement range	800mm ~ 1400mm	
Power Supply	AC110V/220V±10%, 50Hz/60Hz	
Coaxial Illumination	≥30000lx	
Power Consumption	12VA	
Fuse	T1.25AL 250V , T2.5AL 125V	
Whole Packing Volume	0.283m³	
Carton Number	2	
Total Weight	61kg	