

SLIT LAMP MICROSCOPE

YZ3

Binocular Handheld Slit Lamp Microscope

The optical path the angular stereo design,the observation is comfortable,the image is clear;

Lighting with LED light source;

Large battery capacity,uniform spot,easy to switch,long service life,patients feel more comfortable;

Total 900g weight,easy to operate and practical enough.



Specifications

Microscope Type	Converging Microscope
Objective Magnification	Two Steps Magnifications
Magnifications	10× , 16×
Diopter Adjustment	±8D
Pupil Distance	50mm ~ 75mm
Working Distance	60mm
Slit Width	0mm ~ 12mm continuously adjustable
Slit Angle	±30°
Diameter of Light Spot	Φ12mm , Φ5mm , Φ1mm , Φ0.2mm
Filters	Color Temperature Compensation, Redfree, Cobalt blue
Light Source	LED
Maximal Illumination	≥30000Lx
Power Consumption	7.4Vd.c Secondary Li-ion Battery
Working Duration	1.5 hours for common use (max illumination)
Packing Volume	275mm×243mm×255mm
Net Weight	900g