

SLIT LAMP MICROSCOPE

YZ5J

Slit Lamp Microscope

Two-step magnifications : 10×, 16×.

All the optical lenses are moisture-proof, mildew-proof and anti-reflected treated.

Excellent optical performances.

Convenient and compact.



Specifications

Microscope Type	Converging Stereoscopic Microscope
Objective Magnification	Two Steps Magnifications
EyePiece	10×
Magnifications (Field of view)	10×(Φ18mm) 16×(Φ11.25mm)
Pupil Distance	55mm ~ 82mm
Diopter Adjustment	±6D
Slit Width	0mm ~ 9mm adjustable (slit is round when the slit width is 9mm)
Slit Height	1mm ~ 8mm adjustable
Slit Angle	0° ~ 180° adjustable
Slit Inclination	5° , 10° , 15° , 20° four steps
Filters	Heat Absorption, Grey, Redfree, Cobalt blue
Diameter of Light Spot	Φ9mm, Φ8mm, Φ5mm, Φ3mm, Φ2mm, Φ1mm, Φ0.2mm
Power Supply	AC110V/220V±10%, 50Hz/60Hz
Packing Volume	670mm×620mm×490mm
Total Weight	23Kg
Optional Accessories	R-type, T-type Applanation Tonometer, Measuring EyePiece, Trihedron Lens, Gonioscope, Funduscope, Aspheric Lens