

SLIT LAMP MICROSCOPE

YZ5X

Slit Lamp Microscope

Advanced converging stereoscopic microscope.
14mm aperture provides broader and clearer field of view.
Compact and flexible structure ensures the convenient use.



Specifications

Microscope Type	Converging Stereoscopic Microscope
Objective Magnification	Two Steps Magnifications
EyePiece	10×
Magnifications (Field of view)	10×(Φ18mm) 16×(Φ11.25mm)
Pupil Distance	55mm ~ 82mm
Diopter Adjustment	±6D
Slit Width	0mm ~ 14mm adjustable(slit is round when the slit width is 14mm)
Slit Height	1mm ~ 14mm adjustable
Slit Angle	0° ~ 180° adjustable
Slit Inclination	5° , 10° , 15° , 20° four steps
Filters	Heat Absorption, Redfree, Cobalt blue
Diameter of Light Spot	Φ14mm , Φ9mm , Φ5.5mm , Φ0.3mm
Power Supply	AC110V/220V±10% , 50Hz/60Hz
Packing Volume	680mm×550mm×475mm
Total Weight	20Kg
Optional Accessories	Applanation Tonometer, Trihedron Lens, Gonioscope, Funduscope, Aspheric Lens, Teaching Tube