

SLIT LAMP MICROSCOPE

YZ5X1

Slit Lamp Microscope

- Galilean Stereoscopic Microscope Provides the practitioners unsurpassed clarity and sharpness, and avoid double images.
- Five-step drum magnification.
- Optional accessories : Teaching Tube, CCD camera and Tonometer.



Specifications

Microscope Type	Galilean Stereoscopic Microscope
Objective Magnification	Five Steps Drum Magnifications
EyePiece	12.5×
Magnifications (Field of view)	6×(Φ33mm) 10×(Φ22.5mm) 16×(Φ14mm) 25×(Φ8.8mm) 40×(Φ5.5mm)
Pupil Distance	55mm ~ 75mm
Diopter Adjustment	-5D ~ +3D
Slit Width	0mm ~ 14mm adjustable (slit is round when the slit width is 14mm)
Slit Height	1mm ~ 14mm adjustable
Slit Angle	0° ~ 180° adjustable
Slit Inclination	5° , 10° , 15° , 20° four steps
Filters	Heat Absorption, Redfree, Cobalt blue
Diameter of Light Spot	Φ14mm , Φ9mm , Φ5.5mm , Φ0.3mm
Power Supply	AC110V/220V±10% , 50Hz/60Hz
Packing Volume	680mm×550mm×475mm
Total Weight	20Kg
Optional Accessories	Applanation Tonometer, Trihedron Lens, Gonioscope, Funduscope, Aspheric Lens, Teaching Tube