

## Accessories for the Slit Lamp Microscope



Trihedron lens



Gonioscope-Single Lens



Tetrahedron Lens



Contact Lens

Germany glass with multiple coatings provides provides clear image and broad vision.

### Specifications

Image Magnification	2.97×
Field of View	50°
Working Distance	47mm
Clear Aperture	48mm



### Specifications

Image Magnification	0.77×
Field of View	74°
Working Distance	10mm
Clear Aperture	26mm



### Specifications

Image Magnification	0.65×
Field of View	74°
Working Distance	10mm
Clear Aperture	18mm

